

**Imagine...** reliable identification  
without manual labelling



SAKURA

**Tissue-Tek<sup>®</sup> AutoWrite<sup>®</sup>**

Cassette / Slide Printer

# Again a product to make your life easy!



## Tissue-Tek® AutoWrite® NEXT Cassette & Slide Printer

With the **new and improved** version of the Tissue-Tek® AutoWrite®, the **Tissue-Tek® AutoWrite® NEXT**, Sakura Finetek innovates the way to print cassettes and glass slides in the histology and cytology laboratory. Start using the **Tissue-Tek® AutoWrite® NEXT** in your laboratory today and experience the convenience of fully automated and reliable labelling of your cassettes and glass slides... just by connecting them to a LIMS or a stand alone computer.

### Tissue-Tek® AutoWrite® offers you:

- Unmatched printing quality with the new and improved printing system
- Different cassettes or glass slides printed simultaneously
- 1- or 2-dimensional bar coding for easy and reliable identification of cassettes or glass slides
- Convenient loading with the stacked cassettes, offered in all colours and types
- Fast unloading with an unload module
- Continuously printing during the day



Unmatched printing quality by patented ink in a convenient cartridge



Easy and reliable identification with 1-D and 2-D bar coding



Stacked cassettes with easy removable tape for convenient loading



Unloading module for 100 cassettes or 110 glass slides



Load up to 480 cassettes or 450 glass slides at once

## Dedicated to Histopathology

Another product that makes the life of a histotechnologist easy. Full flexibility, high throughput, high quality results, reduction of errors, all in one instrument, offered by Sakura Finetek. As the innovative company in histopathology, Sakura Finetek is continuously looking for possibilities to improve the laboratory... and succeeds in offering a solution for the problems found in the histopathology laboratory.

## Product Specifications

<b>Throughput</b>	Up to 900 cassettes per hour Up to 800 slides per hour
<b>Dimensions</b> W x D x H	46 x 65 x 94 cm
<b>Weight</b>	29 kg
<b>Capacity</b>	
Cassette Printer	6 magazines, 80 cassettes each
Glass Slide Printer	3 magazines, 150 glass slides each
Unload module	100 cassettes or 110 glass slides
Ink cartridge	60,000 cassettes or glass slides or 3 months (whatever comes first)
<b>Printer Driver</b>	Windows® based
<b>Print Resolution</b>	Up to 360 dpi
<b>Characters</b>	1-D and 2-D barcodes, graphics, alphanumeric and roman characters

Item Code	Product Description
8030	AutoWrite® NEXT Cassette Printer 100-240V
8031	AutoWrite® NEXT Glass Slide Printer 100-240V
8032	AutoWrite® NEXT Unload Station Cassette
8033	AutoWrite® NEXT Unload Station Slide

**First we understand.**  
**Then we innovate.**



**Tissue-Tek® Distributor:**

**Sakura Finetek Europe B.V.**  
The Netherlands  
Phone: +31 (0)71 589 83 00  
Fax: +31 (0)71 589 84 88  
Sakura@sakura.eu  
www.sakura.eu